

Reaktionstester

```
1  OpenWindow(0,200,200,800,600,"Reaktionstester",#PB_Window_SystemMenu|
2                                     #PB_Window_ScreenCentered)
3  ButtonGadget(1,5,0,500,50,"Knopf 1")
4  ButtonGadget(2,5,150,500,50,"Knopf 2")
5  StringGadget(3,5,75,500,50,"")
6  StringGadget(4,5,225,500,50,"")
7
8  If StartDrawing(WindowOutput(0))
9      Box(0,0,800,600,RGB(255,255,255))
10
11  Repeat
12
13      Select WindowEvent( )
14
15          Case #PB_Event_CloseWindow
16              End
17
18          Case #PB_Event_Gadget
19              Select EventGadget( )
20                  Case 1
21                      radius=50
22                      Circle(600,170,radius,RGB(255,215,0))
23                      SetGadgetText(3,"Test1")
24
25                  Case 2
26                      SetGadgetText(4,"Test2")
27                      a= ElapsedMilliseconds()
28                      Circle(600,290,radius,RGB(0,255,2))
29              EndSelect
30
31      EndSelect
32
33  ForEver
34      StopDrawing( )
35  EndIf
```

Cursorposition

```
1   Width=800
2   Height=600
3
4   OpenWindow(0,200,200,Width,Height,"Cursorposition",#PB_Window_SystemMenu|#PB_W
5   indow_ScreenCentered)
6       StringGadget(1,5,75,500,50,"")
7
8
9   If StartDrawing(WindowOutput(0))
10      Box(0,0,Width,Height,RGB(255,255,255))
11
12      Repeat
13
14          Select WindowEvent( )
15
16              Case#PB_Event_CloseWindow
17                  End
18
19
20              Default
21                  x=WindowMouseX(0)
22                  y=WindowMouseY(0)
23                  SetGadgetText(1,Str(x)+" | "+Str(y))
24
25          EndSelect
26
27      ForEver
28          StopDrawing( )
29  EndIf
```